Total Body Densitometry – save as *ID_visit#_Bone*




Image not for diagnosis

Printed: 10/07/2020 10:58:53 AM (12.20)76:0.15:153.85:31.2 0.00:-1.00 4.80x13.00 12.8:%Fat=36.6% 0.00:0.00 0.00:0.00 Flename: xhicypc6v.dfb Scan Mode: Standard 0.4 µGy

1 -Statistically 68% of repeat scans fall within 1SD (± 0.010 g/cm² for Total Body Total)

^{2 -}USA (Combined NHANES (ages 20-30) / Lunar (ages 20-40)) Total Body Reference Population

⁽v111) 3 -Matched for Age, Weight (females 25-100 kg), Ethnic

Total Body Densitometry Ancillary Results – save as ID_visit#_BoneExt

Patient: Facility ID:
Birth Date: Referring Physician:
Height / Weight: Measured:
Sex / Ethnic: Analyzed:

ANCILLARY RESULTS [Total Body]

BMD	_				вмс	Area
(g/cm²)	(%)	T-Score	(%)	Z-Score	(g)	(cm ²)
2.435	-	-	-	-	521	214
0.764	-	-	-	-	136	178
1.165	-	-	-	-	434	372
0.847	-	-	-	-	345	408
1.058	-	-	-	-	1,090	1,030
0.771	-	-	-	-	138	179
1.175	-	-	-	-	448	381
0.814	-	-	-	-	355	436
1.131	-	-	-	-	1,287	1,138
0.767	-	-	-	-	274	357
1.170	-	-	_	-	881	753
0.830	-	-	-	-	700	844
0.601	-	-	-	-	199	331
0.983	-	-	-	-	268	273
0.971	-	_	-	-	233	240
1.096	97	-0.4	105	0.7	2,377	2,168
	2.435 0.764 1.165 0.847 1.058 0.771 1.175 0.814 1.131 0.767 1.170 0.830 0.601 0.983 0.971	BMD (g/cm²) (%) 2.435 - 0.764 - 1.165 - 0.847 - 1.058 - 0.771 - 1.175 - 0.814 - 1.131 - 0.767 - 1.170 - 0.830 - 0.601 - 0.983 - 0.971 -	BMD (g/cm²) (%) T-Score 2.435 0.764 1.165 0.847 1.058 0.771 1.175 0.814 1.131 0.767 1.170 0.830 0.601 0.983 0.971	BMD (g/cm²) (%) T-score (%) 2.435	BMD (g/cm²) Young-Adult (%) Age-Matched (%) Z-Score 2.435 - - - - - 0.764 -	BMD (g/cm²) Young-Adult (%) Age-Matched (%) BMC (g) 2.435 - - - 521 0.764 - - - 136 1.165 - - - 434 0.847 - - - 345 1.058 - - - 1,090 0.771 - - - 138 1.175 - - - 448 0.814 - - - 355 1.131 - - - 274 1.170 - - - 881 0.830 - - - - 199 0.983 - - - 268 0.971 - - - 233

^{1 -}Statistically 68% of repeat scans fall within 1SD (± 0.010 g/cm² for Total Body Total)

^{2 -}USA (Combined NHANES (ages 20-30) / Lunar (ages 20-40)) Total Body Reference Population (v111)

^{3 -}Matched for Age, Weight (females 25-100 kg), Ethnic Filename: xhlcypc6v.dfb

Total Body Composition Enhanced Analysis – save as ID_visit#_BCExt

Patient:	Facility ID:
Birth Date:	Referring Physician:
Height / Weight:	Measured:
Sex / Ethnic:	Analyzed:

BODY COMPOSITION

ODT COMP	1		_ 1	_ 1	. 1		1
Region	Tissue (%Fat)	Region (%Fat)	Tissue (g)	Fat (g)	Lean (g)	BMC (g)	Total Mass (kg)
Left Arm	33.7	32.2	2,907	979	1,928	136	3.0
Left Leg	40.3	38.8	11,380	4,589	6,790	434	11.8
Left Trunk	36.8	35.9	14,088	5,182	8,906	345	14.4
Left Total	37.0	35.7	29,745	11,007	18,738	1,090	30.8
Right Arm	33.7	32.1	2,897	975	1,921	138	3.0
Right Leg	40.3	38.8	11,615	4,683	6,933	448	12.1
Right Trunk	36.8	35.9	14,331	5,268	9,063	355	14.7
Right Total	36.2	34.8	31,547	11,428	20,119	1,287	32.8
Arms	33.7	32.2	5,804	1,955	3,849	274	6.1
Legs	40.3	38.8	22,995	9,272	13,723	881	23.9
Trunk	36.8	35.9	28,420	10,450	17,969	700	29.1
Android	42.6	42.1	4,256	1,815	2,441	50	4.3
Gynoid	45.5	44.6	10,754	4,893	5,861	209	11.0
Total	36.6	35.2	61,292	22,435	38,857	2,377	63.7

FAT MASS RATIOS

Trunk/	Legs/	(Arms+Legs)/
Total	Total	Trunk
0.47	0.41	1.07

^{1 -}Statistically 68% of repeat scans fall within 150 (± 0.8 % Fat, ±210 g Tissue Mass, ±520 g Fat Mass, ±610 g Lean Mass for Total Body Total)

Lumbar Spine Ancillary – save as *ID_visit#_LSExt*

Patient: Facility ID:
Birth Date: Referring Physician:
Height / Weight: Measured:
Sex / Ethnic: Analyzed:

ANCILLARY RESULTS [AP Spine]

	BMD 1	Your	Young-Adult Age-Matched 3		вмс	Area	Width	Height	
Region	(g/cm²)	(%)	T-Score	(%)	Z-Score	(g)	(cm²)	(cm)	(cm)
ш	0.993	88	-1.1	91	-0.8	12.03	12.11	3.7	3.26
L2	1.083	90	-1.0	94	-0.6	14.14	13.05	3.8	3.47
L3	1.226	102	0.2	106	0.6	18.76	15.30	4.0	3.80
L4	1.231	103	0.3	106	0.6	20.66	16.78	4.6	3.65
L1-L2	1.040	89	-1.0	93	-0.7	26.17	25.16	3.7	6.72
L1-L3	1.110	95	-0.5	99	-0.1	44.92	40.46	3.8	10.52
L1-L4	1.146	97	-0.3	101	0.1	65.59	57.24	4.0	14.18
L2-L3	1.160	97	-0.3	100	0.0	32.89	28.34	3.9	7.27
L2-L4	1.187	99	-0.1	103	0.3	53.55	45.13	4.1	10.92
L3-L4	1,229	102	0.2	106	0.6	39.42	32.08	4.3	7.46

¹ -Statistically 68% of repeat scans fall within 1SD ($\pm~0.010~\text{g/cm}^2$ for AP Spine L2-L4)

^{2 -}USA (Combined NHANES (ages 20-30) / Lunar (ages 20-40)) AP Spine Reference Population (v102)

^{3 -}Matched for Age, Weight (females 25-100 kg), Ethnic Filename: 8g6kej7r2.dfs

Dual Hip – save as *ID_visit#_DHExt*

Patient: Birth Date: Height / Weight: Sex / Ethnic:

Facility ID: Referring Physician: Measured:

Analyzed:

ANCILLARY RESULTS [DualFemur]

Region	BMD 1 (g/cm²)	Your	g-Adult T-Score	Age-I	Matched Z-Score	BMC (g)	Area (cm²)
Neck Left	1.347		-	-	-	6.47	4.81
Neck Right	1.232	-	-	-	-	5.75	4.67
Neck Mean	1.289	-	-	-	-	6.11	4.74
Neck Diff.	0.115	-	-	-	-	0.72	0.13
Wards Left	1.243	-	-	-	-	3.19	2.57
Wards Right	1.184	-	-	-	-	2.87	2.42
Wards Mean	1.214	-	-	-	-	3.03	2.49
Wards Diff.	0.059	-	-	-	-	0.32	0.14
Troch Left	1.109	-	-	-	-	15.56	14.02
Troch Right	1.055	-	-	-	-	14.22	13.48
Troch Mean	1.082	-	-	-	-	14.89	13.75
Troch Diff.	0.055	-	-	-	-	1.34	0.54
Shaft Left	1.497	-	-	-	-	20.32	13.57
Shaft Right	1.382	-	-	-	-	19.70	14.25
Shaft Mean	1.440	-	-	-	-	20.01	13.91
Shaft Diff.	0.115	-	-	-	-	0.62	0.68
Total Left	1.307	-	-	-	-	42.35	32.40
Total Right	1.224	-	-	-	-	39.68	32.40
Total Mean	1.266	-	-	-	-	41.01	32.40
Total Diff.	0.083	-	-	-	_	2.68	0.00

^{1 -} Statistically 68% of repeat scans fall within 1SD (± 0.010 g/cm² for DualFemur Total) Filename: pacina_gyqsqn7r2.dfe

R or L Radius – save as *ID_visit#_Rad*

* if you scan both radii you will only get data from a single scan, place the other side scans in a separate input folder and repeat the process and simply cut and past the data into your main data document.

Patient: Facility ID:
Birth Date: Referring Physician:
Height / Weight: Measured:
Sex / Ethnic: Analyzed:

ANCILLARY RESULTS [Right Forearm]

Region	BMD (g/cm²)	Youn	g-Adult T-Score	Age-l (%)	Matched Z-Score	BMC (g)	Area (cm²)
Radius UD	0.296	63	-3.8	73	-2.4	1.01	3.42
Ulna UD	0.237	-	-	-	-	0.61	2.59
Radius 33%	0.897	101	0.1	117	1.5	2.07	2.31
Ulna 33%	0.798	-	-	-	-	1.62	2.03
Both UD	0.270	-	-	-	-	1.63	6.01
Both 33%	0.851	-	-	-	-	3.69	4.34
Radius Total	0.613	90	-1.1	102	0.2	8.07	13.17
Ulna Total	0.472	-	-	-	-	5.18	10.96
Both Total	0.549	-	-	-	-	13.25	24.13

¹ -Statistically 68% of repeat scans fall within 1SD ($\pm~0.020~\text{g/cm}^2$ for Right Forearm Radius 33%)

^{2 -}USA (Combined NHANES (ages 20-30) / Lunar (ages 20-40)) Forearm Reference Population (v111)

^{3 -}Matched for Age, Ethnic 9 -Lunar calibration in use.

Filename: palmev_fqigw07r2.dfa